



**KSA POST-OP SHOE**

Made from a soft brush nylon upper and has 8mm micro sole unit. The shoe is supplied with two pieces of hook fastening. Available in three sizes.

Code	Size	Sole Length
FPOSTS	Small	23cm
FPOSTM	Medium	27cm
FPOSTL	Large	30cm



**DARCO® HEEL WEDGE**

Protects heel by significantly reducing weight-bearing pressure. Fits left or right foot.

Code	Size	Men's	Women's
RA63075XS	X-Small	-	up to 4
RA63075S	Small	up to 5	4-5½
RA63075M	Medium	5-7½	6-8
RA63075L	Large	8-10	8½+
RA63075XL	X Large	11+	-

**DARCO® ORTHOWEDGE™**

Protects the forefoot by removing pressure from the metatarsal head. Fits left or right foot.

Code	Size	Men's	Women's
RA6090	X-Small	-	up to 4
RA6091	Small	up to 5	4-5½
RA6092	Medium	5-7½	6-8
RA6093	Large	8-10	8½+
RA6094	X Large	11+	-



**DARCO® APB™ HI-TOP SHOE**

The APB™ for all season closed toe protection. Room for bulky, post-op dressings or even most fibreglass casts. Breathable upper provides high-top protection. Ankle strap holds shoe snug. Opens flat for foot insertion. Fits left or right.

Code	Size	Men's	Women's
RA60327XS	X-Small	-	up to 4½
RA60327S	Small	4-5½	5-7½
RA60327M	Medium	6-7½	8+
RA60327L	Large	8-10	-
RA60327XL	X Large	11+	-



**DARCO® SHOE**

It's soft, padded mesh upper with open toe is lightweight and breathable. Adjustable closures accommodate bulky dressings and provide compression. Fits left or right foot.

Code	Size
RA6931	Small Men: 5-7
RA6932	Medium Men: 7½-9
RA6933	Large Men: 9½-11
RA6934	X Large Men: 11+
RA6951	Small Women: 4
RA6952	Medium Women: 4½-6
RA6953	Large Women: 6½-8



**CANVAS ROCKER CAST SHOE**

Popular style canvas duck upper. Rocker sole helps promote a natural gait while the open toe and heel allows for more versatility in sizing. Elastic strap around the ankle reduces slippage. Navy Blue.

Code	Size	Sole Length
RA64433	Navy Blue, Child, X-Small	7¼" (190mm)
RA64434	Navy Blue, Small	9½" (242mm)
RA64435	Navy Blue, Medium	11¼" (285mm)
RA64436	Navy Blue, Large	12¾" (323mm)
RA64437	Navy Blue, X-Large	13½" (341mm)
RA64438	Navy Blue, Small/Medium	10½" (267mm)
RA64439	Navy Blue, Medium/Large	12¼" (310mm)
RA64440	Teddy Bears	5¼" (133mm)
RA64441	Teddy Bears	7½" (190mm)



**GAITKEEPER™ CAST SHOE**

With cushioned air sole design. Prevents excessive cast wear and reduces impact loading of the injured leg. Enhanced rocker bottom permits natural gait while flexible sole delivers superior comfort. Elastic straps hold cast shoe firmly in position.

Code	Size	Mens up to 4½	Womens up to 5
RA64735	Extra Small	Mens 5-7	Womens 5½-7½
RA64736	Small	Mens 7½-9	Womens 8-9½
RA64737	Medium	Mens 9½-11	Womens 10-11½
RA64738	Large	Mens 11½+	Womens 12+
RA64739	Extra Large		



**DARCO® SLIMLINE CAST BOOT (Adult and Child)**

Has a stiff sole. May be used as a post-op. Long wearing polyester and EVA upper. Rocker design sole with a ripple finish. Velcro® closures. Blue. Sold individually.

Adult Code	Size	Men's	Women's
RA64716S	Small	up to 4	up to 6
RA64716M	Medium	5-7	6½-8
RA64716L	Large	7½-10	8½+

Paediatric Code	Size	Sole Length
M40DRCBPS	Small	14cm
M40DRCBPM	Medium	17cm
M40DRCBPL	Large	20cm



**DARCO® MEDICAL SURGICAL SHOE™**

Strapless forefoot closure eliminates buckle pressure and can expand to accommodate dressings up to 17" in circumference. Ankle strap and raised heel help counter heel slippage, a common complaint with post op footwear. Includes a rocker sole and removable EVA insole that can be modified to eliminate pressure under the surgery site. *Fits right or left.*

Code	Size	Description
RA64457	Black, Womens	Small up to 4
RA64458	Black, Womens	Medium 4½-6
RA64459	Black, Womens	Large 6½-8
RA64460	Black, Mens	Small 5-7
RA64461	Black, Mens	Medium 7½-9
RA64462	Black, Mens	Large 10½-11
RA64463	Black, Medium	Mens X-Large 11½-14
RA64464	Black, Youth	1-3½



**DYNA WALK™**

The DynaWalk™ orthotic system uses adjustable metatarsal bars that can be precisely positioned to offload either the forefoot or the rearfoot. The rigid footplate on the DynaWalk provides a solid platform to significantly reduce flexion/extension and other motion to the foot. Ample size and smooth insert help reduce friction.

Code	Size	Men's	Women's
RA64129XS	Extra Small	-	up to 4
RA64129S	Small	-	4½-6
RA64129M	Medium	5-7	6½-8
RA64129L	Large	7½-9	8½+
RA64129XL	X-Large	9½-11	-
RA64129XXL	XX-Large	11½+	-



**ROYCE® PRESSURE RELIEF SHOE**

For the healing foot. Resilient insole has a series of removable hexagons for tailored pressure relief. Lightweight. Velcro® closure. Rocker bottom. Removable toe piece. Fits left or right.


Code	Size	Men's	Women's
RA61168	Small	up to 5	4½-6½
RA61169	Medium	5½-7½	7-9
RA61171	Large	8-10	9+
RA62865	Extra Large	10½+	-



**DM SYSTEMS CAST WALKER®**

Is the solution for walking cast wearers who are tired of heavy, bulky, costly ambulatory aids and cast protection. CastWalker increases patient mobility and comfort, making a walking cast less of an impediment to an active lifestyle. Patient's gait pattern improves because CastWalker allows for a more normal swingthrough motion. Textured grid on bottom minimizes slipping.

Code	Size
RA64497	Small 6" (15.2cm)
RA64498	Medium 7" (17.7cm)
RA64499	Large 8" (20.3cm)



**MULTI PODUS SYSTEM®**

Multi Podus® System helps correct footdrop contracture and deformity when used in conjunction with an effective ROM program. Applies 10 to 25 lbs. of continuous counterforce pressure to plantar surface of foot. Helps prevent pressure sores by eliminating pressure on heel and keeping bed linens off toes. Rotator bar corrects internal / external rotation of hip, foot and ankle. Detachable outsole provided for transfer from floor to bed. Toe cap included. Lightweight Kydex®. Washable Kodel® liner.

Code	Product	Size
RA62381	Multi Podus System - S	6" - 7½" (15.2cm-19cm)
RA62382	Multi Podus System - M	7½" - 9" (19cm-22.8cm)
RA62383	Multi Podus System - L	9" - 10" (22.8cm-25.4cm)
RA62384	Multi Podus System - XL	10" - 11" (25.5cm-28cm)



**DARCO® WOUND CARE SHOE SYSTEM**

Has a unique insole layering system recommended for the treatment of foot problems where it is desirable to redistribute weight away from specific areas.

Code	Size	Men's	Women's
RA63217S	Small	5-6½	4-6½
RA63217S/M	Small/Medium	7-7½	7-7½
RA63217M	Medium	8-8½	8-8½
RA63217L	Large	9-10½	9+
RA63217XL	X-Large	11-12½	-
RA63217XXL	XX-Large	13+	-



**SEMI RIGID POST-OP SHOE**

Features all the support and comfort of higher priced models. Fits left or right.

Code	Size	Men's	Women's
RA64161	X-Small	-	up to 4
RA64163	Small	up to 5½	4-5½
RA64165	Medium	6-8	6-8½
RA64167	Large	8½-11	9+
RA64168	X Large	11½+	-



**DARCO® WALKER FX PRO™**

Darco® Walker FX Pro™ Ultra low profile design eliminates height difference between walker and normal shoe without altering gait. Padded liner is vented promoting air circulation, keeping the lower leg cool and comfortable. Anatomically contoured "spatula" stirrups comfortably immobilize the ankle while the pressure-moldable, conforming insole absorbs shock. Available in standard and low heights.

Code	Product	Mens	Womens
RA64580	Darco FX Pro Low XS	up to 4	up to 4
RA64581	Darco FX Pro Low Sm	4½-6½	4½-7
RA64582	Darco FX Pro Low Med	7-9½	7½-9½
RA64583	Darco FX Pro Low Lg	10-12	10+
RA64584	Darco FX Pro Low XL	13+	-
RA64585	Darco FX Pro Standard XS	up to 4	up to 4
RA64586	Darco FX Pro Standard Sm	4½-6½	4½-7
RA64587	Darco FX Pro Standard Med	7-9½	7½-9½
RA64588	Darco FX Pro Standard Lg	10-12	10+
RA64589	Darco FX Pro Standard XL	13+	-

**PROSTEP™**

One-piece design. Economical choice for serial casting.

Code	Size	Shoe Size
RA62997S	Small	Men up to 5, Women 5
RA62997M	Medium	Men 5-9, Women 5-9
RA62997L	Large	Men 9+, Women 9+



**SAMSON® SHORT LEG WALKER**

A cost-effective alternative to casting that provides comfort and support. Hyperblended foam liner wraps around leg for exceptional comfort. A natural, crutch-free gait is permitted by the ultra low heel and rocker bottom. The foot is stabilised by unique criss-cross straps at the forefoot. Stays provide excellent support and attach with hook and loop, allowing a customised, reproducible fit. Fits left or right foot.

Code	Size	Men	Women
RA60781	X-Small	5	up to 3
RA60782	Small	5-7	3-5
RA60783	Medium	7-10	5-8
RA60784	Large	10-12	8-10
RA60785	X Large	12+	10+

**EQUALISER® AIR WALKER**

Pneumatic system increases stability. Air-fill level is easily adjusted to achieve desired compression and custom fit. Accommodates the swelling patterns that occur throughout rehab.

Code	Size	Men's	Women's
RA64103XS	Extra Small	up to 3	up to 4
RA64103S	Small	3½-6	4-6
RA64103M	Medium	6½-9½	6½-9½
RA64103L	Large	9½-11½	9½-11½
RA64103XL	Extra Large	11½+	11½+



### ROYCE® DIABETIC WALKER

With insole of removable, resilient hexagons for tailored pressure relief. Fully padded inner bootie protects the foot. Pre-inflated bladders for strut protection. Rocker bottom. Plantar pressure relief. Alternative to total-contact casting 90° ankle fixation reduces shearing

Code	Size	Men's	Women's
RA60137S	Small	up to 7	up to 8
RA60137M	Medium	7½-9½	8½-10
RA60137L	Large	10+	11+



### ACHILLES TENDON WALKER

Instantly immobilises the foot and ankle. Tri-level insole is constructed of Evazote® on the top layer, a 5cm hind foot wedge and a NickelPlast™ base. Adjustable hind foot wedge initially positions the foot at 30° of plantar flexion (50mm high, this may vary with foot length). As recovery progresses, decrease plantar flexion in 2.5mm increments by simply peeling a layer away from the base of the wedge. Progressive adjustments are possible in seconds.

Code	Size
RA63175	Medium Walker, up to size 9
RA63176	Large Walker, size 9 and above
RA63282	Replacement Insole Medium
RA63283	Replacement Insole large



### WALKER LINING

Designed to save the expense of buying new ankle walking braces. This lining replaces older ones and just slips into place. Fastened together by hook and loop.

Code	Size
FLINING	One Size Fits All




### CAM WALKER

Fix the ankle or prescribe a range-of-motion. Features adjustable, double-action joints, fixed or up to 45° ROM. Lock-in dorsiflexion or plantar flexion. Small: 9¾"L x 3¾"W (24.7cm x 9.5cm), Medium: 10¾"L x 4¼"W (27.3cm x 10.7cm), Large: 11½"L x 4½"W (29.2cm x 11.4cm). Replacement liners available

Code	Size	Mens	Womens	Code	Size
RA6048	Small	4-6	4-6	RA60454	Liner Small
RA62445	Medium	6½-9½	6½-9½	RA62446	Liner Medium
RA6049	Large	10+	10+	RA60455	Liner Large

# THOMAS BED SPLINTS

Used to extend the leg in the treatment of fractures. Available in adult and child sizes. For special sizes please state internal ring circumference and inside rod length. Chrome plated. Please specify size when ordering



Code	Description
RTHOMASSR	Split ring, all sizes
RTHOMASFR	Full ring, all sizes
RTHOMASFR REF	Full ring refurbished, all sizes
RTHOMASFR WL	Full ring with loop, all sizes
RTHOMASR REF	Split ring refurbished, all sizes
RTHOMAS RWL	Split ring with loop, all sizes

SUGGESTED USAGE CHART	
Small Child	16" (41cm) Ring, 30" (76cm) Length
Child	20" (51cm) Ring, 34" (86cm) Length
Youth	23" (58cm) Ring, 39" (99cm) Length
Adult	26" (71cm) Ring, 42" (107cm) Length
Large Adult	28" (71cm) Ring, 42" (107cm) Length

Note: Suggested sizes displayed. Any size available (maximum length 42")

