



KSA^{LTD}

Childrens
Product Brochure

Welcome

to our new catalogue

We pride ourselves on being one of the most competitively priced providers of physiotherapy and rehabilitation products. Thanks to our associations with companies such as Markell, Alimed & Benik, our already comprehensive range of products has now even greater depth.

Our full range of products can be accessed via our website: www.ketteringsurgical.co.uk

When ordering please quote:

Official order number
Invoice and delivery Address
Product Code(s) and requirements
Telephone, fax and e-mail details in case of query.

Post: Kettering Surgical Appliances Ltd
73 Overstone Road
Northampton
NN1 3JW

Telephone: 01604 622886

Fax: 01604 629689

E-mail: orders@ketteringsurgical.co.uk

Returns policy:

Made to measure products are non-returnable.
We are happy to accept goods for return or exchange providing they are clean, unused, in their original packaging and in a re-saleable condition.
Goods returned for exchange or credit may be subject to a re-stocking fee.

Payment:

Prepayment is required by non account customers and for private orders.
Payment can be made by cheque or bacs.
Cheques should be made payable to Kettering Surgical Appliances Ltd.

Childrens Products *index*

| | |
|-------------------------------|-------|
| MARKELL CHILDRENS FOOTWEAR | 3-6 |
| WRIST SPLINTS FOR CHILDREN | 7-13 |
| ELBOW ORTHOTICS | 15 |
| KNEE ORTHOSIS | 16 |
| TRUNK ORTHOSIS | 17 |
| VELCRO & HOOK & LOOP FASTENER | 18 |
| SPLINTING MATERIALS | 19-21 |
| INDEX | 22 |



TARSO MEDIUS 1644/2644

A symmetrical straight last boot made with soft tan leathers, a well defined heel pocket, soft rolled top line and special built in padding at the back of the boot and under the tongue. These boots are made in a full range of sizes, from 0000 for newborns, up to size 9 for the older 3-4 year child. Available in narrow or medium width.

| Code | Size |
|--------|---------------------------------------|
| TM1644 | 0000, 000, 00, 0, 1 – 6 in half sizes |
| TM2644 | 6.5 - 9 in half sizes |

Open toe boots have two plated steel screw receptacles set flush in the bottom of each shoe so that splints or appliances can be fastened directly to the shoe bottom without rivets or clamps.



TARSO MEDIUS 1645

The new model 1645 is a Tarso Medius Straight Last Boot, which is both open at the toe, and has a cut-out for the heel. This newly designed boot is made with the younger clubfoot patient in mind. Young infants with atypical or hyper sensitive feet, can sometimes experience difficulty in adjusting to the boots and bars after serial casting. For these young clubfoot patients, Markell have designed the 1645 boot with a special cut out in the heel, to ensure that the heel is down in the correct position, to help relieve pressure, and also to better hold the foot and prevent slipping, pulling out, and blistering. Model 1645 is made in sizes 0000 to 1, and in 3 different widths Narrow, Medium and Wide, to accommodate the full range of infants with slender, normal, or robust feet.

| Code | Sizes |
|--------|---------------------------|
| TM1645 | 0000, 000, 00, 0 + 1 Only |




TARSO PRONATOR

Made on a special full abduction last, which has a more pronounced curvature than our closed toe boots. The full abduction last is made to provide maximum correction for the patient with Metatarsus Adductus, and also for the corrected clubfoot. Available in narrow or medium width.

| Code | Size |
|--------|-------------------------------------|
| TP1934 | 0000, 000, 00, 0, 1-5 in half sizes |

For sizing please see pages 5 and 6




SUEDE HEEL GRIPS

Used to obtain finer fit around the top of the upper of the DB Boots. Sold in pairs.

| Code | Size |
|-----------|---------------------------|
| RHEELGRIP | One Size Fits All (Pairs) |

BOOTS FOR CLUB FOOT


Tarso Open Toe Boots are designed to be used to help maintain correction in babies with Clubfoot and Metatarsus Adductus. Tarso Open Toe Boots are made on both our Tarso Medius straight last and the Tarso Pronator full abduction last.



STANDARD SPLINT M140 DB

Dennis Browne Type – With Foot Plates. Made of quality 24ST aluminium alloy. Gold anodised bar. Bar sizes 6",7",8",9",10" and 12". All foot plates have two large pre-drilled holes for splint adaptor screws and have protractors ruled in 5° increments and clearly labeled with the most widely used settings of 45° and 70°.

| Code | Size |
|------------|--------------|
| M140DB 6" | 6" (15.2cm) |
| M140DB 7" | 7" (17.7cm) |
| M140DB 8" | 8" (20.3cm) |
| M140DB 9" | 9" (22.8cm) |
| M140DB 10" | 10" (25.4cm) |
| M140DB 12" | 12" (30.4cm) |



TARSO ADJUSTABLE BAR SPLINT WITH PLATES 240 AJP

Plates and 1 inch knurled nuts are the same as for 140-DB Standard splint.

Medium: 8"-12" (20.3cm - 30.4cm)
 Large: 12"-16" (30.4cm - 40.6cm)

| Code | Size |
|----------|-----------------------------------|
| M240AJPM | Medium - 8"-12" (20.3cm - 30.4cm) |
| M240AJPL | Large - 12"-16" (30.4cm - 40.6cm) |



AJP ADJUSTABLE BAR 265AJP

The length of the bar can be adjusted by using the allen key supplied. The nut can be tightened with the spanner included.

Medium: 7"-11" (17.7cm - 27.9cm)
 Large: 9"-15" (22.8cm - 38.1cm)

| Code | Size |
|-------------|-----------------------------------|
| M265AJP7-11 | Medium - 7"-11" (17.7cm - 27.9cm) |
| M265AJP9-15 | Large - 9"-15" (22.8cm - 38.1cm) |




BENDING IRONS

Used for shaping the ends of bars. 140DB, 240AJP, 265AJP. Purchase in pairs for best results.

| Code | Size |
|-------|-------------------|
| RBEND | One Size Fits All |

MARKELL


SCALE FOR OPEN TOE TARSO SHOES



Directions: Place heel of foot on this line. Determine size line closest to big toe.
For proper shoe size add growth room (usually 1/4 inch)

MARKELL

SCALE FOR OPEN TOE TARSO SHOES



Directions: Place heel of foot on this line. Determine size line closest to big toe.
For proper shoe size add growth room (usually 1/4 inch)



KSA CHILDRENS EVESWELL WRIST BRACE

The Eveswell Wrist Brace has been designed with the help of Paediatric Occupational Therapists. Design features include a wide area between the thumb and forefinger, flaps built into the body of the brace to protect the back of the hand from the fastening straps plus a soft malleable palmer bar.

Note: KSA Eveswell wrist braces are also available in adult sizes

| Size | Left | Right | Wrist Circ. |
|---------|----------|----------|-------------|
| B SMALL | WCEVEBSL | WCEVEBSR | 8.5cm |
| X SMALL | WCEVEXSL | WCEVEXSR | 10cm |
| SMALL | WCEVESL | WCEVESR | 12cm |
| MEDIUM | WCEVEML | WCEVEMR | 13cm |
| LARGE | WCEVELL | WCEVELR | 14cm |
| X LARGE | WCEVEXLL | WCEVEXLR | 15cm |



KSA CHILDRENS EVESWELL ULNA DEVIATION WRIST BRACE

The Eveswell Ulna Deviation wrist brace has the same features as the Eveswell Wrist Brace but also has a malleable ulna bar.

Note: KSA Eveswell Ulna Deviation wrist braces are also available in adult sizes.

| Size | Left | Right | Wrist Circ. |
|---------|----------|----------|-------------|
| B SMALL | WCEUDBSL | WCEUDBSR | 8.5cm |
| X SMALL | WCEUDXSL | WCEUDXSR | 10cm |
| SMALL | WCEUDSL | WCEUDSR | 12cm |
| MEDIUM | WCEUDML | WCEUDMR | 13cm |
| LARGE | WCEUDLL | WCEUDLR | 14cm |
| X LARGE | WCEUDXLL | WCEUDXLR | 15cm |



KSA CHILDRENS THUMB INHIBITION SPLINT

The Thumb Inhibition Splint helps to push the thumb away from the palm by means of a neoprene pad, which pushes on the thenar eminence.

| Size | Left | Right | Wrist Circ. |
|---------|----------|----------|-------------|
| B SMALL | WCTISBSL | WCTISBSR | 8.5cm |
| X SMALL | WCTISXSL | WCTISXSR | 10cm |
| SMALL | WCTISL | WCTISR | 12cm |
| MEDIUM | WCTISML | WCTISMR | 13cm |
| LARGE | WCTISLL | WCTISLR | 14cm |
| X LARGE | WCTISXLL | WCTISXLR | 15cm |



KSA CHILDRENS THUMB SPICA

The Thumb Spica fits on the inside of the hand it has a malleable strip shaped to push and control the base of the thumb.

| Size | Left | Right | Wrist Circ. |
|---------|----------|----------|-------------|
| B SMALL | WCSPIBSL | WCSPIBSR | 8.5cm |
| X SMALL | WCSPIXSL | WCSPIXSR | 10cm |
| SMALL | WCSPISL | WCSPISR | 12cm |
| MEDIUM | WCSPIML | WCSPIMR | 13cm |
| LARGE | WCSPILL | WCSPILR | 14cm |
| X LARGE | WCSPIXLL | WCSPIXLR | 15cm |



KSA CHILDRENS WEIGHTED HAND SPLINT

The Weighted Hand Splint is designed to fit on the back of the hand and has an open ended pocket to enable weights to be added to control tremor. The weights are available in 1, 2, or 3oz at extra cost.

| Size | Left | Right | Wrist Circ. |
|---------|----------|----------|-------------|
| X SMALL | WCWEIXSL | WCWEIXSR | 10cm |
| SMALL | WCWEISL | WCWEISR | 12cm |
| MEDIUM | WCWEIML | WCWEIMR | 13cm |
| LARGE | WCWEILL | WCWEILR | 14cm |
| X LARGE | WCWIXLL | WCWEIXLR | 15cm |

N.B. Weights 1, 2 & 3oz are available. Code – WCWEI 1oz, 2oz, 3oz



KSA CHILDRENS MUSGRAVE THUMB SPLINT

The Musgrave Thumb Splint is a thumb abduction splint cut in one piece of neoprene, glued together and then heat taped for strength. The splint can be trimmed at the thumb or the wrist for improved fit.

| Size | Left | Right | Wrist Circ. |
|---------|----------|----------|-------------|
| X SMALL | WCMUSXSL | WCMUSXSR | 10cm |
| SMALL | WCMUSSL | WCMUSSR | 12cm |
| MEDIUM | WCMUSML | WCMUSMR | 13cm |
| LARGE | WCMUSLL | WCMUSLR | 14cm |
| X LARGE | WCMUSXLL | WCMUSXLR | 15cm |



CHILDRENS TONE SPLINT

The Neoprene Tone Splint provides low level stretch in pronation or supination. Hook and loop fastening opens down the back of the splint for ease of fitting. The design abducts the thumb, and the strap, which fastens further up the arm, is attached on the back of the hand for easier palm control.

| Size | Left | Right | Wrist Circ. |
|---------|-----------|-----------|-------------|
| X SMALL | WCTONEXSL | WCTONEXSR | 10cm |
| SMALL | WCTONESL | WCTONESR | 12cm |
| MEDIUM | WCTONEML | WCTONEMR | 13cm |
| LARGE | WCTONELL | WCTONELR | 14cm |
| X LARGE | WCTONEXLL | WCTONEXLR | 15cm |



CHILDRENS WRIST BRACE WITH THUMB

This Wrist Brace with Thumb is designed for children where the therapist is unsure whether to fit a wrist brace or a thumb brace. Immobilises the thumb and wrist and helps keep the thumb in the 'penguin' position.

| Size | Left | Right | Wrist Circ. |
|---------|---------|---------|-------------|
| X SMALL | WCNTXSL | WCNTXSR | 10cm |
| SMALL | WCNTSL | WCNTSR | 12cm |
| MEDIUM | WCNTML | WCNTMR | 13cm |
| LARGE | WCNTLL | WCNTLR | 14cm |
| X LARGE | WCNTXLL | WCNTXLR | 15cm |



HAND SIZING KIT

The kit consists of eight graduated sizes of Neoprene gloves with thumb (thumb abduction splints). The splints in this kit are product BD-88 with hook and loop closure. They allow you to physically try on the splint to determine correct size. Please use the sizing kit for any Benik Paediatric Hand products you wish to order.

| |
|----------|
| Code |
| RBPHSKVC |

BD-88 NEOPRENE GLOVE WITH THUMB

Constructed of 3mm neoprene and terry lined for comfort. This glove is designed to offer support for the thumb CMC and MCP. If more support is needed, see thermoplastic options on page 10. Modifications available (see below left). Indicate left or right hand.

| Code | Size |
|-------|--------------|
| BD-88 | AAAA, AAA, |
| | AAB, AA, |
| | A ,B, C or D |

Note: Sizes AAAA and AAA have to be velcro™ closed



BD-88 MOD, OPTION A (VOLAR VIEW)

The glove features the addition of proximal length, palmer pocket with stay and a full length closure.

A variety of options can be added to the BD-88 and RG-87 to address the specific needs of the patient. You can add: PROXIMAL and/or DISTAL length (additional length from sizing kit)

- Dorsal Pocket
- Palmer Pocket
- Aluminium Stay
- Full Length Hook/Loop Closure
- (Partial Closure also available with added proximal length)
- Thermoplast Stay



RG-87 NEOPRENE GLOVE WITH THUMB & INDEX FINGER

Constructed of 3mm neoprene, terry lined for comfort. Like the BD-88 this glove offers support to the thumb CMC, MCP while including the index for support and isolation. Modifications available (See left).

Indicate left or right hand. See thermoplastic options on page 10.

| Code | Size |
|-------|--------------|
| RG-87 | AAAA, AAA, |
| | AAB, AA, |
| | A, B, C or D |



STANDARD THUMB STAY - VOLAR VIEW
Gives extra support for thumb abduction.



STANDARD THUMB STAY - DORSAL VIEW
Gives extra support for thumb abduction.



THUMB WEB SPACE
Web space opener



DORSAL THUMB STAY
Control MCP hyperextension



ULNA STAY
Control Ulna deviation



THENAR SUPPORT
Gives support to the CMC while leaving the MCP and IP free.

THERMOPLASTIC STAYS

Thermoplastic stays can be added to most of Benik's splints and garments where additional support is desired. By warming the thermoplastic stay in a microwave, the stay becomes pliable and contours to your patient's shape and build. The above illustrations (with the thermoplast featured in yellow) show some of the most requested locations for the thermoplastic stays. The chart below indicates thermoplast options available for hand splints.

To order something other than pictured, please provide a tracing of the hand and outline the desired coverage area on both dorsal and volar sides.

| THERMOPLAST AVAILABILITY | BD-88 | BD-88 MOD OPTION A | RG-87 | RG-87 MOD OPTION A | W-303 | W-302 |
|--------------------------|-------|--------------------|-------|--------------------|-------|-------|
| Standard Thumb Stay | ✓ | ✓ | ✓ | ✓ | ✓ | N/A |
| Thumb Web Space | ✓ | ✓ | ✓ | ✓ | ✓ | N/A |
| Dorsal Thumb Stay | ✓ | ✓ | ✓ | ✓ | ✓ | N/A |
| Ulnar Stay | N/A | ✓ | N/A | ✓ | ✓ | ✓ |
| Thenar Support | ✓ | ✓ | ✓ | ✓ | ✓ | N/A |



W-302 PAEDIATRIC WRIST SUPPORT WITH STAYS

Constructed of 3mm Neoprene and terry lined for comfort, this wrap-around design gives maximum support for extension by incorporating dorsal and volar pockets with removeable aluminium stays. Secured by a wrap-around D-ring wrist strap and hook/loop closures located on the distal and proximal end. (see thermoplastic options available on page 10)

| Code | Size |
|-------|----------------------------------|
| W-302 | XXS, XS, SM, MED, LG or XL |



W-303 WRIST SUPPORT WITH THUMB

Constructed with the same features as the W-302, the W-303 provides additional neoprene support for the CMC and MCP joints of the thumb. (See thermoplastic options available on page 10.)

| Code | Size |
|-------|----------------------------------|
| W-303 | XXS, XS, SM, MED, LG or XL |

W-300 Wrist Series

Refer to chart below and select Paediatric size XXS, XS, SM, MED, LG or XL. (Indicate left or right hand)

| STOCK PAEDIATRIC SIZE | SIZING KIT SIZE | SPLINT LENGTH MCP's to FOREARM | APPROXIMATE AGE |
|-----------------------|-----------------|--|-----------------|
| XXS | AAAA | 3 ³ / ₄ " (9.5cm) | <1 |
| XS | AAA | 4 ¹ / ₄ " (10.8cm) | <2 |
| SM | AAB/AA | 5 ¹ / ₄ " (13.5cm) | 2-3 |
| MED | A | 6" (15.25cm) | 4-5 |
| LG | B/C | 6" (15.25cm) | 6+ |
| XL | D | 7" (17.5cm) | 10+ |

For custom sizes, refer to Hand Splint Order Form on page 13.



W-211 WRIST WRAP

This neoprene wrist wrap has a 1/16" (1.5mm) mouldable thermoplastic stay that completely surrounds the CMC and MCP joints to aid in positioning and immobilising the thumb.

| Code | Size |
|-------|--|
| W-211 | AAAA, AAA, AAB, AA, A, B, C or D |



W-312 WRIST SUPPORT

Thermoplastic version of Benik's W-302 wrist support (page 11). A 1/16" (1.5mm) thermoplast panel replaces dorsal and volar stays which help control Ulna deviation, flexion and hyperextension of the wrist. See stock sizing chart on page 11.



W-313 WRIST SUPPORT WITH NEOPRENE THUMB

Thermoplastic version of Benick's W-303 wrist support. The dorsal and volar stays have been replaced with a 1/16" (1.5mm) mouldable thermoplast panel to help control Ulna deviation, flexion and hyperextension of the wrist. The Neoprene extends to the IP joint of the thumb while the stay extends over the CMC to the MCP joint of the thumb. See stock sizing chart on page 11.



W-323 WRIST SUPPORT WITH NEOPRENE AND THERMOPLAST THUMB

An adaptation to the W-313 by extending the thermoplast towards the IP joint of the thumb. See stock sizing chart on page 11.



PEDI COMFY™ HAND/WRIST SPLINT

For support and progressive extension of the paediatric hand and wrist. Foam-covered metal skeleton bends to fit a wide range of problems. Fits right or left hand.

| Code | MP's | Wrist Circumference |
|-----------|---------------------|----------------------|
| RA510431S | 0"-3" (0cm-7.5cm) | 2½" (6.5cm) |
| RA510431M | 2½"-3½" (6.5cm-9cm) | 2½"-3" (6.5cm-7.5cm) |
| RA510431L | 4-5" (10cm-12.5cm) | 3½"-4½" (9cm-11.5cm) |



PEDI COMFY™ FUNCTIONAL HAND THUMB ORTHOSIS

Removeable finger separators can be positioned to prevent pressure areas between PIP and DIP joints in any angle. Fits right or left hand

| Code | MP's | Wrist Circumference |
|-----------|---------------------|----------------------|
| RA510433S | 0"-3" (0cm-7.5cm) | 2½" (6.5cm) |
| RA510433M | 2½"-3½" (6.5cm-9cm) | 2½"-3" (6.5cm-7.5cm) |
| RA510433L | 4-5" (10cm-12.5cm) | 3½"-4½" (9cm-11.5cm) |



PUCCI® R.I.P. PEDI WHFO

Maintains functional position. Removable, heat-mouldable plastic platform supports fingers in neutral and inhibits spasms of finger joints. Includes closed-cell foam splint, foam palmar roll, platform, one cover and mesh laundry bag. One size. Specify Left or Right.

| Code | Size |
|----------|-------|
| RA510861 | Left |
| RA510861 | Right |



PUCCI® EZE PEDI HAND ORTHOSIS

Incorporates long opponens and elongated "C" bar for maintenance of first web space and normal arches of hand. Maintains functional wrist position. Bendable to address wrist contractures. Includes closed-cell foam splint, foam palmar roll, one splint cover and mesh laundry bag. One size. Specify Left or Right.

| Code | Size |
|----------|-------|
| RA510859 | Left |
| RA510859 | Right |



ELBOW WRAP WITH PAD

Elbow pad with adjustable closures. Contoured hook and loop sensitive elbow wrap with an extra layer of 3mm neoprene covering the posterior elbow. Adjustable closures at forearm and bicep allow for maximum adjustment.

| Code | Lateral Length | Accommodates Elbow Circ. |
|-----------|----------------|--------------------------|
| E400VPed3 | 3" (76mm) | 5"-7" (127mm-177mm) |
| E400VPed4 | 4" (102mm) | 6"-8" (152mm-203mm) |
| E400VPed5 | 5" (127mm) | 7"-8" (177mm-203mm) |
| E400VPed6 | 6" (152mm) | 8"-9" (203mm-228mm) |
| E400VPed7 | 7" (177mm) | 9"-11" (228mm-280mm) |



HINGED ELBOW WRAP

The dual pivot aluminum hinge (with removable 165° extension stop) in this versatile neoprene wrap provides light medial/lateral support and stability to the elbow. Adjustable dual closure design allows for maximum custom fit with a stock size.

Options include removable 1"(25mm) elastic strapping that can be added to offer additional resistance to hyperextension - if needed specify desirable length. Proximal and distal clinch straps are removeable.

| Code | Lateral Length | Circumference |
|----------------|----------------------------|----------------------|
| E505VPed6 | 6" (152mm) | 5"-10" (127mm-254mm) |
| E505VPed7 | 7" (177mm) | 5"-10" (127mm-254mm) |
| E505VPed8 | 8" (203mm) | 6"-10" (152mm-254mm) |
| E505VPed9 | 9" (228mm) | 6"-11" (152mm-280mm) |
| PELS El Straps | Optional additional straps | |



PEDI COMFY™ ELBOW SPLINT

A padded maleable wire splint that can be bent to many positions to fit a wide range of sizes. Very useful for progressive splinting. Comes with two sets of covers and laundry bag. Fits either arm.

| Code | Size |
|---------|-------------------------------|
| 510432S | Small fits - average 1-4 yrs |
| 510432M | Medium fits - average 4-8 yrs |
| 510432L | Large - 9 yrs to small adult |

| Elbow Splint Size | Bicep Circ. | Forearm Circ. | Length |
|-------------------|----------------------------|----------------------------|-------------|
| Sm (purple) | 0"- 6"(<15cm) | 0"- 6"(<15cm) | 5½"(14cm) |
| Med (turquoise) | 6½"-7" (16.5cm-17.5cm) | 6½"-7" (16.5cm-17.5cm) | 6½"(16.5cm) |
| Lg (green) | 7"-8" (17.5cm - 20.5cm) | 7"-8" (17.5cm - 20.5cm) | 7"(16.5cm) |



PUCCI® INFLATABLE PEDI ELBOW

Treats flexion contractures with constant low stretch of contracted tissue. Air-stay allows some movement, preventing unforgiving pressure against soft tissue. For severe contractures, an elbow pad is available. Sizing: Measure circ. 3" (7.5cm) above and below elbow. If upper arm is 4"-6½" (10-16.5cm) and forearm is 3½"-5½" (9-14cm) order small 6"L (15cm). If measurements are larger, measure 4" (10cm) above and below elbow joint. If upper arm is 7"-9" (17.8-22.8cm) and forearm is 5"-7" (12.7-17.8cm) order regular 8"L (20.3cm).

| Code | Description |
|------------|--|
| RA510864S | Upper arm 4"-6½" (10-16.5cm), Forearm 3½"-5½" (9-14cm) |
| RA5100864L | Upper arm 7"-9" (17.8-22.8cm), Forearm 5"-7" (12.7-17.8cm) |
| RA510865 | Elbow Pad |



PAEDIATRIC KNEE WRAP WITH PAD

Contoured Hook and loop sensitive knee wrap with an extra 3mm layer of neoprene covering the anterior surface. Wrap design allows for maximum adjustment and future growth. Open popliteal for comfort. Custom sizes also available.

| Code and Size | Lateral Length | Mid-Patella Circ. |
|---------------|----------------|-----------------------|
| K400V Ped 3 | 3" (76mm) | 5" – 8" (127-203mm) |
| K400V Ped 4 | 4" (102mm) | 7" – 10" (177-254mm) |
| K400 V Ped 5 | 5" (127mm) | 9" – 12" (228-305mm) |
| K400 V Ped 6 | 6" (152mm) | 10" – 14" (254-355mm) |
| K400 V Ped 7 | 7" (178mm) | 11" – 15" (280-380mm) |



PAEDIATRIC HINGED KNEE WRAP

The dual pivot aluminium hinge (with removable 165° extension stop) in this versatile neoprene hinged paediatric knee wrap provides light medial and lateral support to ligaments and offers stability to weak knees. Open popliteal and patella offer comfort and reduction of bulk. Options include removable 1" (25mm) elastic straps that can be added to offer additional resistance to hyper-extension if needed. Specify desired length. Custom sizes also available.

| Code and Size | Lateral Length | Mid-Patella Circ. | Distal Circ. |
|---------------|----------------|-------------------|-------------------|
| K305V Ped 6 | 6" (152mm) | 5-11" (127-280mm) | 5-9" (127-228mm) |
| K305V Ped 8 | 8" (203mm) | 6-14" (152-355mm) | 5-11" (127-280mm) |
| K305V Ped 10 | 10" (254mm) | 7-15" (177-380mm) | 6-12" (152-305mm) |
| K305V Ped 12 | 12" (305mm) | 8-16" (203-406mm) | 6-14" (152-355mm) |

Pels El. Straps



PAEDIATRIC BASIC KNEE SLEEVE

The Basic Knee Sleeve is constructed of 3mm terrycloth lined neoprene. This sleeve offers compression and warmth. Custom orders only. See below for measurements required. Specify colour.

| Code - K-100T | Measurements Required |
|---------------|---|
| G | Mid Patella Circumference |
| HD | Using Anterior Side – Mid patella to distal end |
| J | Circumference at distal end |
| HP | Using Anterior Side – Mid patella to proximal end |
| I | Circumference at proximal end |



PEDI COMFY™ KNEE SPLINT

A foam-encased malleable wire splint that bends to accommodate any angle and fit a wide range of sizes. Very useful for progressive splinting.

| Code | Knee Splint Size | W&L | Strap Circ. |
|----------|------------------|--------------------|-------------------|
| RA61061S | Sm (purple) | 5"x6¼" (12.5x16cm) | 0"-16" (0-40.5cm) |
| RA61061M | Med (turquoise) | 5¼"x7" (13.5x18cm) | 0"-16¼" (0-41cm) |
| RA61061L | Lg (green) | 6¼"x8¼" (16x21cm) | 0"-20" (0-51cm) |